AT-2

6-Channel Electrocardiograph
6-Kanal-Elektrokardiograph
Electrocardiographie à 6 canaux

THE ART OF DIAGNOSTICS

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**Article Number 2. 510199e
Sept.95, March 96, June 2000,
May 2001, August 2001**

Associated Documents

Guide to the SCHILLER Interpretation and Measurement Program E/ D/ F
Article Number 2. 510 179

SCHILLERAG

Altgasse 68

6341 Baar, Switzerland

Phone: ++ 41 41 766 42 42

Fax: ++ 41 41 761 08 80

www.schiller.ch

CE 0123

93/42/EEC Medical Devices:
0123 'Notified Body' TÜV P. S.

Where to Obtain Service and Sales Advice
Kundendienst und Verkaufs-/Beratungsstellen
Coordonnées de nos services d'assistance technique et commerciale

USA / Canada	WELCH ALLYN SCHILLER Inc. , 7420 Carroll Road, San Diego, CA , US-92121-2334 USA	
	Tel.: +1 858 635 6023	Fax : +1 858 635 6611
	Home Page : www.welchallyn.com	
Asia Pacific	SCHILLER Asia Pacific , 10 Jalan SS 3/33, Taman Universiti, 47300 Petaling Jaya, Selangor, Malaysia	
	Tel.: + 603 7877 5336	Fax : + 603 7877 5744
Austria	SCHILLER HmbH , Kampmüllerweg 24, A-4044 Linz, Austria	
	Tel.: + 43 732 709 90	Fax : + 43 732 757 000
France	SCHILLER Medical S.A , BP 50, 19, Avenue de la Gare, F-67162 Wissembourg / Cedex, France	
	Tel.: +33 3 88 63 36 00	Fax : +33 3 88 94 12 82
Germany <i>(EU authorized representative)</i>	SCHILLER Medizintechnik GmbH , Rudolf-Diesel Strasse 14, D-85521 Ottobrunn, Germany	
	Tel.: + 4989 629 981 0	Fax : + 4989 609 509 0
India	SCHILLER Healthcare India Pvt. Ltd. , D.C.Silk Mills Compound, 'A' Wing, 1st floor, 5, Chunawala Estate, Kondivitta Lane, Andheri - Kurla Road, Andheri (E, Mumbai - 400 059, India	
	Tel.: + 9122 826 3520	Fax : + 9122 826 3525
Italy	ESAOTE Spa (SCHILLER) , Via di Caciolle 15, I-50127 Firenze, Italy	
	Tel.: + 39055 4229 201	Fax : + 39055 4229 208
Switzerland	SCHILLER Reomed AG , Riedstrasse 14, CH-8953 Dietikon, Switzerland	
	Tel.: +411 744 3000	Fax : + 411 740 3710
All other countries	SCHILLER AG , Altgasse 68, CH-6341 Baar, Switzerland	
	Tel.: + 4141 766 4242	Fax : + 4141 761 0880
	Home Page : www.schiller.ch	

This equipment has been tested and found to comply with the limits for a class A digital device, pursuant to both Part 15 of the FCC (Federal Communications Commission) Rules and the radio interference regulations of the Canadian Department of Communications. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with this instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Disposal Instructions and Battery Care

Battery Care and Disposal Instructions

- DO NOT DISPOSE OF THE BATTERY BY FIRE OR INCINERATOR - DANGER OF EXPLOSION
- DO NOT ATTEMPT TO RECHARGE THE BATTERY - DANGER OF EXPLOSION
- DO NOT OPEN THE BATTERY CASING - DANGER OF ACID BURN

*Only dispose of the battery in official recycling centres or municipally approved areas.
Alternatively, used batteries can be returned to Schiller AG for disposal.*

Unit Disposal Instructions

*Units no longer required can be returned to Schiller AG for disposal. Alternatively
dispose of the unit in municipally approved recycling centres.*

Power Supply

The mains connection is on the rear of the unit.

The power supply voltage is set by the factory for 100-115V (nom. 110V) or 220-240V (nom. 230V) working. The setting is indicated by the indented metal strip on the fuse panel. Contact your dealer if the voltage needs to be changed.

The mains indicator lamp on the keyboard is always lit when the unit is connected to the mains supply. The unit can either be operated from the mains supply or from the built-in rechargeable battery.

Changing a Mains Fuse

If it is necessary to change a fuse, always replace with the correct rating i.e 2x200mAT for 230V, or 2x315mAT for 110V .

To change a fuse press the two retaining lugs on side of the fuse panel (situated below the mains connector on the back panel. Remove the fuse panel and replace the fuse(s). Click back the fuse panel.

Intended Use

The CARDIOVIT AT-2 is a 6-channel ECG recorder with all (12) ECG signals simultaneously processed to provide instant ECG recordings. Two automatic recording modes can be individually preset to enable one button ECG recording of preferred print formats.

Terms of Warranty

The SCHILLER AT-2 is warranted against defects in material and manufacture for the duration of one year (as from date of purchase). Excluded from this guarantee is damage caused by an accident or as a result of improper handling. The warranty entitles free replacement of the defective part. Any liability for subsequent damage is excluded. The warranty is void if unauthorized or unqualified persons attempt to make repairs.

In case of a defect, contact your dealer or the manufacturer.

The manufacturer can only be held responsible for the safety, reliability, and performance of the apparatus if:

- assembly operations, extensions, readjustments, modifications, or repairs are carried out by persons authorized by him, and
- the AT-2 and approved attached equipment are used in accordance with the manufacturers instructions.

THERE ARE NO EXPRESS OR IMPLIED WARRANTIES WHICH EXTEND BEYOND THE WARRANTIES HEREINABOVE SET FORTH. SCHILLER MAKES NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WITH RESPECT TO THE PRODUCT OR PARTS THEREOF.

About this Handbook

The philosophy of SCHILLER is one of continuous improvement. Our aim is to provide the user with the most up-to-date information and the latest technological developments.

Your suggestions and comments are welcome on all SCHILLER documentation. Please contact the SCHILLER Technical Documentation Department.

PHYSICIAN'S RESPONSIBILITY

THE AT-2 ELECTROCARDIOGRAPH IS PROVIDED FOR THE EXCLUSIVE USE OF QUALIFIED PHYSICIANS OR PERSONNEL UNDER THEIR DIRECT SUPERVISION. THE NUMERICAL AND GRAPHICAL RESULTS FROM A RECORDING MUST BE EXAMINED WITH RESPECT TO THE PATIENTS OVERALL CLINICAL CONDITION. THE RECORDING PREPARATION QUALITY AND THE GENERAL RECORDED DATA QUALITY, WHICH COULD EFFECT THE REPORT DATA ACCURACY, MUST ALSO BE TAKEN INTO ACCOUNT.

IT IS THE PHYSICIANS RESPONSIBILITY TO MAKE THE DIAGNOSIS OR TO OBTAIN EXPERT OPINION ON THE RESULTS, AND TO INSTITUTE CORRECT TREATMENT IF INDICATED.

FEDERAL LAW IN THE USA RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

Safety Notices

TO PREVENT ELECTRIC SHOCK DO NOT DISASSEMBLE THE UNIT. NO SERVICEABLE PARTS INSIDE. REFER SERVICING TO QUALIFIED PERSONNEL ONLY.

DO NOT USE THIS UNIT IN AREAS WHERE THERE IS ANY DANGER OF EXPLOSION OR THE PRESENCE OF FLAMMABLE GASES SUCH AS ANAESTHETIC AGENTS.

IN THE EVENT OF ACCIDENTAL LCD BREAKAGE AND RESULTANT LEAKAGE OF FLUID, DO NOT INHALE, INGEST OR MAKE CONTACT WITH THE SKIN. IF CONTACT IS MADE RINSE IMMEDIATELY.

THIS PRODUCT IS NOT DESIGNED FOR STERILE USE.

THIS PRODUCT IS NOT DESIGNED FOR OUTDOOR USE.

SWITCH THE UNIT OFF BEFORE CLEANING AND DISCONNECT FROM THE MAINS.

DO NOT, UNDER ANY CIRCUMSTANCES, IMMERSE THE UNIT OR CABLE ASSEMBLIES IN LIQUID.

THE DEVICE MUST ONLY BE OPERATED USING BATTERY POWER IF THE EARTH CONNECTION IS SUSPECT OR IF THE MAINS LEAD IS DAMAGED OR SUSPECTED OF BEING DAMAGED.

DO NOT USE HIGH TEMPERATURE STERILISATION PROCESSES (SUCH AS AUTOCLAVING). DO NOT USE E-BEAM OR GAMMA RADIATION STERILISATION.

DO NOT USE SOLVENT CLEANERS

USE ONLY ACCESSORIES AND OTHER PARTS RECOMMENDED OR SUPPLIED BY SCHILLER AG. USE OF OTHER THAN RECOMMENDED OR SUPPLIED PARTS MAY RESULT IN INJURY INACCURATE INFORMATION AND/ OR DAMAGE TO THE UNIT.

THE AT-2 COMPLIES WITH EMC REGULATIONS FOR MEDICAL PRODUCTS WHICH AFFORDS PROTECTION AGAINST EMISSIONS AND ELECTRICAL INTERFERENCE. HOWEVER SPECIAL CARE MUST BE EXERCISED WHEN THE UNIT IS USED WITH HIGH FREQUENCY EQUIPMENT.

IT MUST BE ENSURED THAT NEITHER THE PATIENT NOR THE ELECTRODES (INCLUDING THE NEUTRAL ELECTRODE) COME INTO CONTACT WITH OTHER PERSONS OR CONDUCTING OBJECTS (EVEN IF THESE ARE EARTHED).

THERE IS NO DANGER WHEN USING THE ECG UNIT FOR A **PACEMAKER PATIENT** OR WITH SIMULTANEOUS USE OF OTHER ELECTRICAL STIMULATION EQUIPMENT. HOWEVER, THE STIMULATION UNITS SHOULD ONLY BE USED AT A SUFFICIENT DISTANCE FROM THE ELECTRODES. IN CASE OF DOUBT, THE PATIENT SHOULD BE DISCONNECTED FROM THE RECORDER.

THIS UNIT IS CF CLASSIFIED ACCORDING TO IEC 601-1. THIS MEANS THAT THE PATIENT CONNECTION IS FULLY ISOLATED AND DEFIBRILLATION PROTECTED. SCHILLER CAN ONLY GUARANTEE PROTECTION AGAINST DEFIBRILLATION VOLTAGE HOWEVER, WHEN THE ORIGINAL SCHILLER PATIENT CABLE IS USED.

IF SEVERAL UNITS ARE COUPLED THERE IS A DANGER OF SUMMATION OF LEAKAGE CURRENT

DO NOT TOUCH THE CASING DURING DEFIBRILLATION

IF THE PATIENT CABLE SHOULD BECOME DEFECTIVE AFTER DEFIBRILLATION, LEAD OFF WILL BE DISPLAYED AND AN ACOUSTIC ALARM GIVEN

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
THE ART OF DIAGNOSTICS

Contents - User Guide AT-2

INTRODUCTION	2	SETTINGS FOR AUTOMATIC MODE	22
OPERATION MODES	4	Average Cycles	23
Automatic Mode	4	Measurements and Markings	24
Manual Mode	5	Interpretation	25
LOCATION & POWER	6	Interpretation Settings	26
Location	6	Selecting Rhythm Leads	27
Power Supply	6	PATIENT CABLE CONNECTIONS	28
BASIC INFORMATION	7	Standard Leads	30
Switching On and Off	7	RECORDING AN ECG IN AUTOMATIC MODE	31
Potential Equalisation	7	RECORDING AN ECG IN MANUAL MODE	33
KEYBOARD	8	CARE & MAINTENANCE	35
INDICATORS	10	Care of your CARDIOVIT AT-2	35
GENERAL SETTINGS	11	Self-test	35
Defaults	12	12 Monthly Check	36
Language	15	Cleaning the Print Head	36
Baseline Filter	16	REPLACING THE RECORDING PAPER	38
Mains Filter	17	TROUBLE SHOOTING	39
Myogram Filter	18	ORDERING INFORMATION	40
Defining Lead Sequence & Printout	19	TECHNICAL DATA	41
Acoustic QRS Indication	20		
Time / Date	21		

Introduction

Introduction

The CARDIOVIT AT-2 is a 6-channel ECG recorder. All ECG signals are simultaneously processed to provide instant ECG recordings. Two automatic recording modes can be individually preset to enable one button ECG recording of preferred print formats .

Individual lamps are provided to give power, paper error, filter, lead group and lead off indications.

In addition, any detected disturbance (i.e. loose electrode or end of paper), gives an audible alarm and the corresponding indicator lamp flashes.

The AT-2 includes the following features:

- Low weight and compact dimensions
- Large A4 size printout from integrated quality thermal printer
- Built-in rechargeable battery for mains-independent use
- Simple one key operation
- Automatic or manual recording modes
- Selectable printing formats
- ECG memory for easy copying
- Interpretation program option (including measurements)

Operation Modes


Automatic Mode

Automatic Mode provides a printout giving 10 seconds of ECG recording of all 12 leads in 2 different formats.

The following can be programmed freely in each of the formats before recording:

- Lead Format
- Chart Speed
- With the optional interpretation program installed it is also possible to select the measurement table, average cycles with optional markings and interpretation statements for the printout.

For further information see paragraph "Settings for Automatic Mode".


Manual Mode

Manual Mode provides a real time printout of 6 leads that are selected and indicated on the keyboard.

The following can be freely selected before or during recording:

- Lead Group
- Chart Speed
- Sensitivity
- Myogram Filter

For further information see paragraph "ECG Recording in Manual Mode".

Location & Power

Location

Do not keep or operate the apparatus in a wet, moist, or dusty environment. Also, avoid exposure to direct sunlight or heat from other sources. Do not allow the unit to come into contact with acidic vapours or liquids, as such contact may cause irreparable damage. The unit should not be placed near X-ray or diathermy units, large transformers or motors.

The unit must be placed on a flat surface. The unit should not be operated in areas where there is any danger of explosion.

Power Supply


The unit can either be operated from the built-in rechargeable battery, or from the mains.

The mains connection is on the rear of the unit. The mains indicator lamp is always lit when the unit is connected to the mains supply.

A battery indicator lamp confirms battery operation. When the battery capacity is limited, the indicator flashes. To recharge the battery, connect the apparatus to the mains supply by means of the supplied power cable. A totally discharged battery needs less than 15 hours to be fully recharged (60% in less than 3 hours, 90% in less than 7 hours).

A fully charged battery lasts approximately 6 hours of normal use. The unit can remain connected to the mains supply without any danger of damage to either the battery or the unit.

Switching On and Off

The CARDIOVIT AT-2 is switched on by means of the green  key and off by means of the red  key.


The unit is switched off after 5 minutes (30 seconds if battery capacity is limited) if no key is pressed and the patient cable is not connected.

Potential Equalisation









If the AT-2 is used in conjunction with other patient connected equipment, we recommend that the potential equalisation stud () on the rear of the unit is connected to the hospital/ building common ground with the yellow/green ground cable (Part-no. 2.310005).

When working from an emergency vehicle, the vehicle common ground can be used.

Keyboard

Keyboard

	Switch unit on		Cancel or enable QRS beeper
	Switch unit off		Automatic ECG sensitivity adjustment from 10 to 5 mm/mV in order to avoid overlapping traces in automatic mode only
	Copy stored ECG		
	Start automatic recording	  	ECG sensitivity selector (5, 10 or 20 mm/mV) manual recording only
	Start manual recording	  	Chart speed selector (5, 25 or 50 mm/s) manual recording only
	Stop recording/move paper to start position		
	Switch myogram filter (muscle tremor filter) on or off		Key for 1 mV indication mark on output during manual recording. Use this key also for baseline recentering
	Lead group selector (forward)		Key for initiation of setups and selection of second format for printout
	Lead group selector (backward)		

Indicators

○ ~ Mains indicator (lit when mains connected)

○  Battery lamp (lit when running on battery power - mains not connected) (blinking when battery capacity is limited)

○ LEAD OFF Warning lamp for loose electrode connection or poor electrode contact

○ PAPER Warning lamp for end of paper or paper jam

○ FILT Myogram filter (lit when filter ON)

I II III
aVR aVL aVF ○ Indicator lamp for selected lead group (Standard)
(Cabrera: aVL, I, -aVR, II, aVF, III) in manual mode only

V1 V2 V3
V4 V5 V6 ○ Indicator lamp for selected lead group (Standard) (in manual mode only)


II aVF III
V2 V4 V5 ○ Indicator lamp for selected lead group in manual mode only

V4 V5 V6
V7 V8 V9 ○ Indicator lamp for selected lead group in manual mode only

General Settings

Each parameter is set by means of a code. This code comprises a combination starting with the **ALT** key followed by a number of keys and is always confirmed with the **STOP** key. As soon as the **ALT** key is pressed, the keyboard is dedicated to the programming function.

The setting is remembered and the keyboard is only released for other functions when the **STOP** key is pressed. Once the settings have been confirmed, they are stored in the memory even when the unit is switched off. As an example, if you want to set the language on your AT-2 to English, proceed as follows:

On the following pages the programmable parameters and the programming sequences are described in detail.

NOTE THAT THE SETTINGS ARE ALWAYS CONFIRMED BY PRESSING THE  KEY.

General Settings

Defaults

Default Settings		
Setting	AT-2	AT-2 with Interpretation Option
Language	English	English
Leads	Standard (S)	Standard (S)
Format 1	ECG: 25mm/s short (o)	ECG: 25 mm/s, short (o)
		MECG: 2 x 6 (25 mm/s) + 2
		Measurements: suppressed (-)
		Marks: enabled (+)
		Interpretation: enabled (+)
Format 2	ECG: 25 mm/s, long (ooo)	ECG: 25 mm/s, long (ooo)
		MECG: none (-)
		Measurements: disabled (-)
		Marks: enabled (+)
		Interpretation: disabled (-)
Rhythm Leads	V1 = ECG; V1, II = MECG	V1 = ECG; V1, II = MECG
Autom. Centering	Enabled (+)	Enabled (+)
Printout of Signals	Sequential	Sequential
Baseline Filter Setting	0.05 Hz	0.05 Hz
Mains Filter Setting	50 Hz (USA - 60Hz)	50 Hz (USA - 60Hz)
Myogram Filter Setting	35 Hz, OFF	35 Hz, OFF
Interpretation Settings		N/A: - U: + A30: - S: -

The defined formats and settings can be checked as follows:

Setup Printout				
Entry Key Sequence				Action
ALT	0	1	1	Printout of programmed settings

A printout of the defined settings will be produced and gives the following information, depending on the installed software:

Unit designation (AT-2),	Software option installed (C = Interpretation) and Software version
Serial number	
Leads	Standard (S) or Cabrera (C)
ECG Format	Long (ooo), Short (o) or Suppressed (-)
MECG	Average cycles as defined in auto ECG recording setup (e.g. 4 * 3 (25 mm/s)+2)
Measurements	Enabled (+) or Suppressed (-)
Marks	Enabled (+) or Suppressed (-)
Interpretation	Enabled (+) or Suppressed (-)
Selected Rhythm leads	Leads selected for R1, R2 resp.
Automatic Centering	Enabled (+) or Suppressed (-)
Printout of signals	Sequential or Simultaneous

General Settings

Baseline Filter	0.05, 0.15 or 0.30 Hz
Mains Filter	50, 60 Hz or OFF (-)
Myogram Filter	25 or 35 Hz, ON (+) or OFF (-)
Interpretation settings:	N/A: +/- ('normal/abnormal' is written (+) or suppressed (-); U: +/- ('unconfirmed report' is written (+) or suppressed (-); A30: +/- (patient age is assumed to be < 30 (-) or >30 (+); S: +/- (low (-) or high (+) sensitivity)

To reset the unit to the basic default settings, proceed as follows:

Reset Setup				
Entry Key Sequence				Action
ALT	0	6	6	Reset to default base settings.

Language

The language is selected as follows:

Select Language					
Entry Key Sequence				Language	Confirm
ALT	0	2	1	German	Press STOP Key
			2	English	
			3	French	
			4	Swedish	
			5	American	
			6	Italian	
			7	Spanish	
			8	Portuguese	

Confirm the selection by pressing **STOP**.

Note: Once selected, the language remains active. However, when the unit is reset to the default settings, the default language will be English.

General Settings

There are three different filters which can be set individually as follows:

- Baseline filter
- Mains filter
- Myogram filter

The setting for each filter is given on the setup printout.

Baseline Filter

The digital **Baseline filter** suppresses excessive baseline drifts. The setting options are as follows:

Baseline Filter				
Entry Key Sequence		Filter Setting		Confirm
ALT	5	0	0.05 Hz (Default)	Press STOP Key
		1	0.15 Hz	
		3	0.30 Hz	

The set value is the lower limit of the frequency range and is normally set to 0.05 Hz. The settings 0.15 and 0.30 Hz should only be used when absolutely necessary, as the possibility exists that they could affect the original ECG signal, especially the ST segments.

Confirm the selection by pressing **STOP**.

Mains Filter

The **Mains filter** is an adaptive digital interference filter designed to suppress AC interference without attenuating or distorting the ECG.

Set the mains filter in accordance with the frequency of your local mains supply as follows:

Mains Filter				
Entry Key Sequence			Filter Setting	Confirm
ALT	8	5	Mains Filter 50 Hz	Press STOP Key
		6	Mains Filter 60 Hz	
		9	Mains Filter off	

General Settings

Myogram Filter

The **Myogram filter** suppresses disturbances caused by strong muscle tremor. The cutoff frequency of the myogram filter is set to 25 or 35 Hz. The myogram filter is switched on and off manually with the **FILT** key. Switching on or off can also be programmed as default when the unit is switched on. When the filter is active, the '**FILTER**' lamp on the unit is lit.

Myogram Filter				
Entry Key Sequence			Setting	Confirm
ALT	8	2	Myogram Filter 25 Hz	Press STOP Key
		3	Myogram Filter 35 Hz	
		1	Myo. Filt. ON when switching on unit (marked on printout with +)	
		8	Myo. Filt. OFF when switching on unit (marked on printout with -)	

Confirm the selection by pressing **STOP**.

An ECG recorded in auto mode is stored unfiltered. It is therefore possible to print the stored ECG either with or without passing the myogram filter.

Defining Lead Sequence & Printout

The required settings can be selected as follows:

Sequences & Print					
Entry Key Sequence			Definition	Confirm	
ALT	7	1	Standard Lead Sequence	Press STOP Key	
		2	Cabrera Lead Sequence		
		3	Simultaneous Print		
		4	Sequential Print		
		5	Auto-Centering ON		
		6	Auto-Centering OFF		

Confirm the selection by pressing **STOP**.

The selectable lead groups for the AT-2 are:

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General Settings

The selectable printout forms are:

Simultaneous	All ECG leads are printed in the same time segment (in automatic mode only).
Sequential	Each group is a contiguous time segment of approximately 2.5 or 5 seconds (in automatic mode only).
Auto-Centering ON	All ECG traces are centred dynamically for optimal use of paper width.
Auto-Centering OFF	ECG traces are set to a fixed baseline position and may possibly overlap.

Acoustic QRS Indication

The acoustic QRS indication can be switched on or off at any time by pressing the key.

Time / Date

CHANGE TIME / DATE						
Time:						
ALT	9	1	1	HHMMSS		beep
Date:						
ALT	9	2	2	DDMMYY		beep
Wintertime to Summertime (+1 Hr)						
ALT	9	4	4			
Summertime to Wintertime (-1 Hr)						
ALT	9	5	5			

Settings for Automatic Mode

Two separate formats for the automatic ECG output can be selected as follows:

Automatic ECG Format		
Entry Key Sequence		Setup Format
ALT	1	Setup Format 1
	2	Setup Format 2

With this the 2 automatic mode formats are defined as detailed on the following pages.

The ECG format is set as follows:

ECG Format					
Entry Key Sequence			Output Format		Confirm
ALT	1 or 2	1	5	No leads are printed	Press STOP Key
			6	Leads are printed in short form (1 sheet)	
			7	Leads are printed in long form (2 sheets)	
			8	Chart speed is 25 mm/s	
			9	Chart speed is 50 mm/s	
			0	Leads are printed in format 4 x 3 + 1 Rhy	

Confirm the selection by pressing **STOP**.

Note: Lead selection for the rhythm lead is defined on page 27.

Settings for Automatic Mode

Average Cycles

The Average cycles are defined as follows:

Average Cycles (Interpretation Option only)				
Entry Key Sequence			Output Format	Confirm
ALT	1 or 2	2	5 No average cycles are printed.	Press STOP Key
			6 4 * 3 (25 mm/s) + 2 * Rhy (25 mm/s) The average complexes are printed out in four groups of three leads with 2 rhythm leads at a chart speed of 25 mm/s.	
			7 4 * 3 (50 mm/s) + 2 * Rhy (25 mm/s) The average complexes are printed out in four groups of three leads with 2 rhythm leads at a chart speed of 50 mm/s.	
			8 2 * 6 (50 mm/s) + 2 * Rhy (25 mm/s) The average complexes are printed out in two groups of six leads with two rhythm leads at a chart speed of 50 mm/s.	

Confirm the selection by pressing **STOP**.

Note: Lead selection for the 2 rhythm leads is defined on page 27.

Settings for Automatic Mode

Measurements and Markings

To define the measurements and markings proceed as follows:

Measurements and Markings (Interpretation Option only)					
Entry Key Sequence			Output Format		Confirm
ALT	1 or 2	3	5	Detailed table of measurement results is omitted (However, the values of electrical axes, intervals, and heart rate are not suppressed.).	Press STOP Key
			6	Detailed table of measurement results is printed.	
			7	Reference markings (beginning and end of P wave and QRS as well as end of T wave) are omitted.	
			8	Reference markings (beginning and end of P wave and QRS as well as end of T wave) are added to ECG cycles.	

Confirm the selection by pressing **STOP**.

Interpretation

To print or suppress interpretation statements on the printout proceed as follows:

Interpretation					
Entry Key Sequence				Output Format	Confirm
ALT	1 or 2	4	5	Interpretation is omitted.	Press STOP Key
			6	Interpretation is printed.	

Confirm the selection by pressing **STOP**.

Full details of the interpretation option are given in the SCHILLER ECG Measurement and Interpretation booklet (art. no. 2.510 179).

Settings for Automatic Mode

Interpretation Settings

The interpretation settings enable the user to determine whether or not certain comments will be added to the interpretation statements on the ECG printout. Furthermore, the patient's age can be defined (< or > 30) and if low or high sensitivity should be applied. Low sensitivity will suppress certain nonspecific ECG diagnosis; this may be advisable when carrying out ECGs for screening.

Interpretation Settings				
Entry Key Sequence			Setting	Confirm
ALT	6	1	'Normal/abnormal' is not printed	Press STOP Key
		2	'Normal/abnormal' is printed.	
		3	'Unconfirmed report' is not printed.	
		4	'Unconfirmed report' is printed.	
		5	Patient age is assumed to be < 30	
		6	Patient age is assumed to be > 30	
		7	Low sensitivity	
		8	High sensitivity	

Settings for Automatic Mode

Selecting Rhythm Leads

The rhythm leads are printed out as defined. Two separate rhythm leads can be selected. The following formats can be set:

Rhythm Leads		
Entry Key Sequence		Setup Lead
ALT	3	Setup Rhythm Lead 1
	4	Setup Rhythm Lead 2

The 2 rhythm leads are defined as follows:

Extremity Leads					
Entry Key Sequence				Lead	Confirm
ALT	3 or 4	8	1	I	Press STOP Key
			2	II	
			3	III	
			4	aVR	
			5	aVL	
			6	aVF	

Precordial Leads					
Entry Key Sequence				Lead	Confirm
ALT	3 or 4	9	1	V1	Press STOP Key
			2	V2	
			3	V3	
			4	V4	
			5	V5	
			6	V6	

Confirm the selection by pressing **STOP**.

Patient Cable Connections

The accessory kit of the electrocardiograph includes a 10-lead patient cable. This cable is plugged into the patient cable socket on the right-hand side of the unit and secured with the two screws.

The CARDIOVIT AT-2 is CF rated. The patient connection is fully isolated and defibrillation protected. Protection against defibrillation voltage is however only ensured, if the original SCHILLER patient cable (Part-no. 2.400070 / USA: 2.400071) is used. Make sure that during ECG recording neither the patient nor the conducting parts of the patient connection or the electrodes (including the neutral electrode) come into contact with other persons or conducting objects (even if these are earthed).


The quality of the ECG is dependent on the preparation and the resistance between the skin and the electrode. To ensure a good quality ECG and minimise the skin/electrode resistance, remember the following points:

1. Ensure that the patient is warm and relaxed.
2. Shave electrode area before cleaning.
3. Thoroughly clean the area with alcohol.
4. Place the **C4** electrode first - in the fifth intercostal space on midclavicular line.
Then place:
 - **C1** in fourth intercostal space at the right sternal border
 - **C2** in fourth intercostal space at the left sternal border
 - **C3** between, and equidistant to, C4 and C2
 - **C6** on left midaxillary line on the same level as C4
 - **C5** between, and equidistant to, C4 and C6

Following these simple guidelines will ensure good results every time.

Patient Cable Connections

Standard Leads

IEC		AHA	
N	Black	RL	Green
R	Red	RA	White
C1	White/Red	V1	Brown/Red
C2	White/Yellow	V2	Brown/Yellow
C3	White/Green	V3	Brown/Green
C4	White/Brown	V4	Brown/Blue
C5	White/Black	V5	Brown/Orange
C6	White/Violet	V6	Brown/Violet
L	Yellow	LA	Black
F	Green	LL	Red

Recording an ECG in Automatic Mode

In **automatic mode**, a full 12-lead ECG is printed in one of two predefined formats with a sensitivity of 10 mm/mV. These two formats are selected by the user to suit his specific needs and requirements (as detailed previously).

When the  key is pressed before recording in automatic mode, the unit detects very large waveform amplitudes and sets the sensitivity for the extremity and/or precordial leads to 5 mm/mV to reduce the overlapping of traces.

To start the automatic ECG recording in Format 1, simply press .

To start the automatic recording in the second format, press  followed by

.

Recording an ECG in Automatic Mode

The printout provides you with the following information:

- ECG recording of all leads in either Standard or Cabrera format according to selection
- Sensitivity
- Heart Rate
- Speed
- Filter Settings
- Patient Data field to manually insert patient data

and if set:

- Average Cycles (Interpretation option only)
- Intervals (Interpretation option only)
- Axis (Interpretation option only)
- Sokolow Index (ECG index for hypertrophy; Interpretation option only)
- Detailed Measurement Table (Interpretation option only)
- Interpretation (Interpretation option only)

- To obtain an extra printout of the ECG recording in Format 1, simply press  .

- To obtain an extra printout of the second format, press  followed by  .

Recording an ECG in Manual Mode

Manual mode provides a direct printout of the real-time ECG with full control of parameter selection.

The following can be freely chosen during or before the recording:

- Lead Group (by means of the  and  keys)

The following lead groups are selectable:

- I, II, III, aVR, aVL, aVF (Cabrera: aVL, I, -aVR, II, aVF, III)
- V1, V2, V3, V4, V5, V6
- II, aVF, III, V2, V4, V5
- V4, V5, V6, V7, V8, V9

- Chart Speed (by means of the ,  and  keys)

- Sensitivity (by means of the ,  and  keys)


- Filter see page 18.

Recording an ECG in Manual Mode

To start the manual recording of a real-time ECG, press the  key.

The printout provides you with the following information:

- The group of the six selected leads with lead identification.
- On the lower edge chart speed, sensitivity and filter settings (if on) are given.
- At the top, the heart rate as current average of 4 beats is shown.

To re-centre the ECG traces, press the  key during operation.

Finish the recording by pressing the  key.

WARNING:

AFTER HEAVY ARTEFACTS OR LEAD OFF, THE INDICATION OF THE HEART RATE
MAY NOT BE RELIABLE.

Care & Maintenance

Care of your CARDIOVIT AT-2

The patient cable should not be exposed to excessive mechanical stress. Whenever disconnecting the leads, hold the plugs and not the cables. Align the leads in such a way as to prevent anyone stumbling over them or any damage caused by the wheels of instrument trolleys. The cable can be wiped with soapy water. Sterilization, if required, should be done with gas only and not with steam. To disinfect, wipe the cable with any standard hospital disinfectant.

The casing of the CARDIOVIT AT-2 should be cleaned with a soft cloth on the surface only.

DISCONNECT THE UNIT BEFORE CLEANING. DO NOT, UNDER ANY CIRCUMSTANCES, IMMERSE THE APPARATUS INTO A CLEANING LIQUID OR STERILIZE WITH HOT WATER, STEAM, OR AIR.

Self-test

Initiate a self-test of the AT-2 as follows:

Initiate Self-Test				
Entry Key Sequence				Action
ALT	0	3	3	Printout of Self-test

A table giving information for the service staff is printed out.

12 Monthly Check

The unit should undergo a technical safety check every 12 months. This safety check should extend to include the following:

- Visual inspection of the unit and cables.
- Electrical safety tests according to IEC 601-1 and IEC 601-2-25.
- Functional tests according to the Service Handbook.

The test results must be documented.

Cleaning the Print Head

If the printer is used a lot, a residue of printers ink (from the grid on the printer paper) can build up on the the print head. This can cause the print quality to deteriorate. We recommend therefore that every month the print head is cleaned with alcohol as follows:

- Remove the paper tray. The printhead is found under, and in from, the paper tray release catch.
- With a tissue dampened with alcohol, gently rub the print head to remove the ink residue. If the printhead is badly soiled, the colour of the grid ink will show on the tissue.

Replacing the Recording Paper


The recording paper must be replaced as soon as the end of the paper is indicated by a red stripe on the lower edge. After the indication first appears, there are about 8 pages left. However, we recommend that the paper be replaced immediately.

If no paper is left, the printing process is interrupted and the paper warning lamp starts to blink. After the paper has been replaced, the printout is restarted by pressing **COPY** or **MAN START** or **AUTO START**.

A step-by-step description of how to change the paper is given on the next page.

Replacing the Recording Paper

- Place fingers under the retaining bar and pull directly upwards. The paper tray cover releases.
- Withdraw the cover from the unit. **DO NOT FORCE, THE PAPER TRAY COVER RUNS FREELY OVER THE DEDICATED RUNNERS.**
- Remove any remaining paper from the paper tray.
- Place a new paper pack into the paper tray with the printed (grid) side facing upwards.
- Place the beginning of the paper over the black paper roller on the paper tray cover.
- Return the paper tray cover in position and press firmly until secure.
- Press the **STOP** key to transport the paper to the start position.

SCHILLER can only guarantee perfect printouts when SCHILLER original chart paper or chart paper of the same quality is used.

Problem	What to Check
Unit does not switch On/ Mains Indicator Lamp is not lit	Check if mains cable is plugged in. Call your local SCHILLER dealer if problem is still present.
"Noisy" traces	CHECK ELECTRODE CONTACT. As much as possible, ensure that patient is relaxed and warm. Activate myogram filter to reduce muscle tremor. Check mains filter to 50 or 60 Hz according to local power supply.
ECG trace 'wanders' away from centre	Baseline drift - check electrode contact. Press 1mV key to reset baseline. Select a higher baseline frequency.
Poor quality printout / ECG traces 'breaking up'	Thermal print head dirty - clean printhead with alcohol. Possible faulty printhead, contact local service centre. Ensure that the paper tray cover is clicked into place
No printout	Connect unit to the mains supply. Ensure that the paper tray cover is clicked into place. Check paper level.

Ordering Information

Your local representative stocks all the disposables and accessories available for the AT-2. In case of difficulty or to obtain the address of your local dealer, please contact the head office. Our staff will be pleased to help process your order or to provide any details for all SCHILLER products.

DESCRIPTION	PART-NO.
10-lead Patient Cable, Standard	2.400070
10-lead Patient Cable, USA	2.400071
Electrodes (box of 500 clip electrodes)	2.155031
Mains (Power) Cable (Germany)	2.300005
Mains (Power) Cable (Switzerland)	2.300003
Mains (Power) Cable (USA)	2.300001
Potential Equalisation (Ground) cable	2.310005
Recording Paper, Z-folded	2. 15 017
User Guide (English / German / French)	2.510 199
User Guide (Italian / Spanish / Portuguese)	2.510 200
Software (C) Interpretation	5. 020002
Guide to the Interpretation and Measurements Programs (English /German / French)	2.510179

Technical data subject to change without notice.

Dimensions	400 x 330 x 100 mm
Weight	4.25 kg (5.05 kg with full paper tray)
Mains Supply	100 to 115 / 220 to 240 VAC, 50/60 Hz
Battery	Built-in 12 V lead-acid battery (rechargeable)
Power Consumption	Recording: 28 VA max
Leads	Standard / Cabrera
Paper Speed	5 / 25 / 50 mm/s (direct)
Sensitivity	5 / 10 / 20 mm/mV, either automatically adjusted or manually selected
Chart Paper	Thermoreactive - Z-folded, 210 mm wide, perforation 280 mm
Printing Process	High-resolution thermal print head, 8 dots per mm
Recording Tracks	6 channels, positioned at optimal width on 200 mm, automatic baseline adjustment
Automatic Lead Programs	6 channel representation of 12 simultaneously acquired standard leads

Technical Data

Data Record: Listing of ECG recording data

Version C: ECG measurement results (intervals, amplitudes, electrical axes), Sokolow Index, average complexes with optional measurement reference markings, and interpretation.

ECG Storage: Memory for 10 s, 12-lead ECG

Circular input memory for 10 s, 12-lead ECG.

Frequency Range of Digital Recorder:

0 to 150 Hz (IEC)

0 to 150 Hz (AHA)

ECG Amplifier:	Simultaneous, synchronous registration of all 9 active electrode signals (= 12 standard leads)
	Sampling frequency: 1000 Hz
	Digital resolution: 5 μ V
	Dynamic range: ± 10 mVAC
	Max. electrode potential: ± 300 mVDC
	Time constant: 3.2 s
	Frequency response: 0.05 to 150 Hz (-3 dB)
	Input impedance: > 10 M Ω
Myogram Filter (muscle tremor filter)	
	25 Hz or 35 Hz, programmable (not active on averaged waveform). The stored ECGs can be printed with or without filter.
Line Frequency Filter:	Distortion-free suppression of superimposed 50 or 60 Hz sinusoidal interferences by means of an adaptive digital filter.
Patient Input:	Fully floating and isolated, defibrillation protected.
Patient Leakage Current:	< 5 μ A

Technical Data

Safety Standard:	CF according to IEC and complying with the following
	RL 93/42/EEC
	EN 60601-1:1990
	IEC 601-1
	IEC 601-2-25:1993
	pr EN 1441:1994
EMC:	CISPR 111: 1985, EN 55011: 1992
	IEC 801-2: 1991
	IEC 801-3: 1984
	IEC 801-4: 1988
	IEC 801-5:
Safety Class:	I according to IEC 601-1 (with internal power supply)
	Ila according to RL 93/42/EEC, CE-0123
Environmental Conditions:	Temperature, Operating: 10° to 40° C
	Temperature, Storage: -10° to 50° C
	Relative humidity: 25 to 95% (non condensing)
	Atmospheric pressure: 700 to 1060 hPa
Control Panel:	Rubber keys

Technical data subject to change without notice.

Available Configurations

The CARDIOVIT AT-2 is available in two versions:

Standard Version:	Unit with ECG recording and printout capabilities.
Version C:	Unit with additional ECG Interpretation program (including measurements).

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